




SMAXEL

## Vaginal treatment mode



Manual LVR Hand piece



Automatic LVR Hand piece

SMAXEL

## Surgical Mode / Dental Mode

### Surgical handpiece



50mm



100mm



E.N.T Hand piece

### Dental Handpiece



SMAXEL

## Fractional Mode



SMAXEL

## Low Pain Fractional Mode



**WORLD GATE**  
Medical & Lab. Supplies




**WORLD GATE**  
Medical & Lab. Supplies





## Fractional (Fractional mode / Low pain fractional mode)



### Various Beam Pattern

- Various beam pattern for each treatment



### Ultra Pulse Beam Minimizes Thermal Damage and Side Effect




### Fractional Surgery indication

Skin tightening	Wrinkle	Large pore	Acne scar
Skin resurfacing	Stretch mark	Laser peeling	Rubbing mode
Pudenda	LVR		



# Smaxel → Surgical Mode / Dental Mode

## Smart mode for Surgical Field (General Sugery , Dental Surgery)



### General Surgery Indication

Corns, Warts etc	Freckles
Birthmarks(mole)	Skin resurfacing
Phimosisctomy	Tumor elimination
Hemorrhoids and anal fistulas	Fine wrinkle
Condyloma	Syringoma
Snoring	

### Dental Surgery Indication

The second implant operation	Hyper sensitive
Hemostasis after dental extraction	Herpetic gingivitis
Crown lengthening	Frenulum excision
Gingival excision, Gingival plastic operation	Cleaning of dental cuticle cyst
Removal of pigmentation	Pulp capping -Fluorine coating
Allergic splenitis	Sterilization in paradental cyst
Sterilization in root canal	Abscess treatment









## **Fractional handpiece + LVR Handpiece**

(Vaginal treatment)


### **Fractional handpiece** (for Labium Treatment)




### **Smart protocol for LVR by Manual Treatment**



### **Fractional handpiece + Auto Vaginal Rejuvenation Handpiece**




## **How** does it work ?



high density laser beam delivers to vagina widely and it makes vaginal tightening & rejuvenation.

## Technical Specification

Model		SMAXEL
Laser Type		CO2 Laser
Wave length		10.600nm
Laser power		1-40W
Function		
Aiming Beam	Diode Laser 655nm (5mw) :1-10 Step	
Surgical Mode	Handpiece Type	50mm,100mm(Basic) E.N.T & Dental Handpiece::Optional
	Laser Operation Mode	1. CW 2. Normal Pulse 3. Super Pulse 4. Super Pulse Continuous 5. Ultra Pulse 6. Smart Protocol(General Surgery/Dental Surgery)
	Pulse Duration	· Ultra Pulse : 90us-900us · Super Pulse :1000us-3500us · Normal Pulse : 4ms-20ms
	Repeat Time	OFF, 2ms/500ms adjustable, 5ms, 10ms, 30ms, 50ms, 70ms, 100ms, 250ms, 500ms
Fractional Mode	Beam Mode	Ultra Pulse
	Density Level	1-23 Step(x,y line) ; 9dots-4489 dots
	Operation Mode	Fractional Mode Low Pain Fractional Mode
	Beam Patterns	Array,Grid, Random
	Beam Size	2x2mm - 20x20mm
	Beam Shape	
	Depth Level	1-5 Step
	Handpiece	· Fractional Handpiece(Galvanic Motor) · Vaginal Rejuvenation Handpiece(Optional)
	Repeat Time	OFF, 0.5 , 1.5, 2, 2.5ms
General Features		
Electrical Requirements		220v, 6A. single phase,50/60hz
Cooling System		water cooling
Dimensions		350x360x1200
Weight		52Kg
Display		10.2 inches color lcd







**SMAVEL**





**WORLD GATE**  
Medical & Lab. Supplies



**Sales; Erbil:** 07503175819

**Dohuk:** 07503032672

**Sulaimanya:** 07508528458

**Service & Maintenance:** 07828386448

**Manager:** 07512354023

[sales@worldgate-medical.com](mailto:sales@worldgate-medical.com)

[www.worldgate-medical.com](http://www.worldgate-medical.com)

**Facebook:** World Gate Medical